## CoFix™ posterior fixation system

### Package contents

The CoFix™ implant consists of one CoFix™ interlaminar implant and two CoFix™ screw/sleeves.

The individual package contains one CoFix™ implant.

### Description

The CoFix™ implant is a spinous process fixation device that stabilizes the spinous processes and spine to act as an adjunct to fusion. The CoFix™ interlaminar implant has two sets of wings. A screw & sleeve are inserted through the holes in the wing and through a prepared hole in the spinous process to fix the crimped wings to the superior and inferior spinous processes.

### Material

The CoFix™ implant is made from wrought titanium 6-aluminum 4-vanadium alloy (ISO 5832-3). It is supplied sterile and is available in numerous sizes.

### Intended use

The CoFix™ implant is a posterior, non-pedicle supplemental fixation device intended for permanent implantation as an adjunct to different fusion techniques in the lumbar spine (L1–S1) in one level. Each wing of the CoFix™ interlaminar implant is attached to the spinous process with the use of a CoFix™ screw/sleeve.

It is intended for attachment to the spinous processes for the purpose of achieving stabilization to promote fusion in the following conditions: Instability, including degenerative disc disease – defined as back pain of discogenic origin with degeneration of the disc, facet joints and pseudospondylolisthesis up to Meyerding grade I, confirmed by history and radiographic studies.

The CoFix™ implant is designed to be used in patients with adequate facet joints as defined by a minimum of 50 % of the medial facet and all facet capsules after the implantation of the CoFix™ implant.

The CoFix™ implant provides additional segmental stabilization to increase the segmental stiffness and promote the fusion process.

The CoFix™ implant is not intended for stand-alone use.

### Indications

It is intended for application at all levels of the lumbar spine (L1–S1) in cases of back pain and discogenic pain due to lumbar degenerative processes of instability as pseudospondylolisthesis (up to Meyerding grade I), osteochondrosis, instability of facet joints and lumbar stenosis or lumbar disc prolapses.

### Contraindications

Do not use the CoFix™ implant in cases of:

- Use in more than one level
- Any medical or surgical condition precluding the potential benefit of spinal surgery
- Spondylolisthesis greater than grade 1 or unstable spondylolisthesis
- Moderate to severe degenerative lumbar scoliosis (Cobb angle of greater than 25°)
- Previous history of osteopenia, osteoporosis, or osteomalacia
- Prior decompressive laminectomy, hemilaminectomy or significant lamina fenestration at the index lumbar level(s) which weakens the spinous process
- Radiographically compromised vertebral bodies at the index lumbar level(s) caused by current or past trauma or tumor (e.g., compression fracture)
- Acute or chronic systemic, spinal or localized infections
- Obesity
- Pregnancy

- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Unwilling or unable to follow postoperative care instructions
- Previous history of allergy to any component of the CoFix™ device
- Previous history of severe allergy or anaphylaxis
- Previous history of disease of bone metabolism
- Previous history of malignancy

### Adverse reactions

Adverse reactions may include:

- Clinical failure (i.e. pain or injury) due to bending, loosening, wear and tear, fracture of implant, loss of fixation, dislocation and/or migration
- Failure of the fusion to heal or spontaneous fusion
- The spine may undergo adverse changes (e.g., adjacent level deterioration) or deterioration including loss of proper spinal curvature, correction, height, and/or reduction, or malalignment, and another surgery may be required
- Pain and/or abnormal sensations due to the presence of the implant
- Primary and/or secondary infections
- Allergic reactions to implant material
- Neurological injury
- Vertebrae and/or spinous process fracture
- Injury to vessels, nerves, dura and organs
- Hematoma and/or impaired wound healing; hemorrhage
- Venous thrombosis, lung embolism, and cardiac arrest
- Death

### Safety precautions

- Prior to use, thoroughly read these instructions for use and become familiar with the surgical technique for implantation of the CoFix™ implant.
- Keep the instructions for use accessible to all staff.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques. Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- The operating surgeon must be capable of performing alternate fusion methods (e.g., pedicle screw fixation) / surgical approaches.
- The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrect operating techniques, the limitations of treatment methods or inadequate asepsis.
- Under no circumstances may modular implant components from different suppliers be combined.
- Each patient's record shall document the implant used (name, article number, lot number).
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed about post-surgical behavioral requirements.
- Damage to the weight-bearing structures can give rise to loosening, dislocation and migration, as well as other complications. To ensure the earliest possible detection of such catalysts of implant dysfunction, the implant must be checked periodically post operative using appropriate techniques.
- Fracture of the CoFix™ implant wing may occur during preparation when the wings are bent open or during insertion when the wings are crimped.

- Never reuse an implant. Although the implant may appear undamaged, previous stresses may have created non-visible damage that could result in implant failure.

- Never use implants if the packaging is damaged. An implant with damaged packaging might be damaged itself and thus may not be used.

- Never use implants that are past their expiration date.

### MRI compatibility

Non-clinical evaluation has demonstrated that the CoFix™ implant is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5-Tesla (1.5T) or 3.0-Tesla (3.0T)
- Spatial gradient field of up to:
  - 7,480 G/cm (74.80 T/m) for 1.5T systems
  - 3,740 G/cm (37.40 T/m) for 3.0T systems.
- Maximum whole body averaged specific absorption rate (SAR) of:
  - 1.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 1.5T
  - 1.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 3.0T.

### 3.0T RF heating

The CoFix™ implant is expected to produce a temperature rise of less than 1.0 °C at a maximum whole body averaged specific absorption rate (SAR) of 1.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 3.0 T Siemens Trio (MRC20587) MR scanner with SYNGO MR A30 4VA30A software as determined by non-clinical evaluation..

### 1.5T RF heating

The CoFix™ implant is expected to produce a temperature rise of less than 3.0 °C at a maximum whole body averaged specific absorption rate (SAR) of 1.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 1.5T Siemens Espree (MRC30732) MR scanner with SYNGO MR B17 software as determined by non-clinical evaluation.

*The RF heating behavior does not scale with static field strength. Devices which do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.*

### MR Artifact

The shape of the image artifact follows the approximate contour of the device and extends radially up to 21 mm from the implant in a 3.0T system with spin-echo sequencing as determined by non-clinical evaluation.

### Storage, inspection and sterilization

#### Storage

The implant is individually packed in protective packaging that is labeled according to its contents.

The implant is sterilized with gamma sterilization (25 kGy minimum).

- Always store the implant in the original protective packaging.
- Do not remove the implant from the packaging until immediately before use.

- Store the implant in a dry and dust-free place (standard hospital environment).

### Disinfection / cleaning

The CoFix™ implant is not designed to be disinfected or cleaned by the user.

### Resterilization

The CoFix™ implant is not designed to be resterilized by the user.

### Procedure

The CoFix™ interlaminar implant must be implanted only with the applicable coflex-F™ trial and instruments. The coflex-F™ trial and instruments are available from the manufacturer at any time. The CoFix™ screw/sleeve must be implanted only with the applicable coflex-F™ instruments. The coflex-F™ instruments are available from the manufacturer at any time.

### Preoperative

The operating surgeon draws up an operation plan specifying and documenting the following:

- Implant component(s) and their dimensions.
- Proper position of the implant component(s) between the spinous processes.
- Determination of intraoperative orientation points.

The following conditions must be fulfilled prior to application:

- All required implant component(s) are readily available.
- Highly aseptic operating conditions are present
- All requisite implantation instruments must be available and in working order, including the coflex-F™ trial and coflex-F™ instruments.

*Never use or process damaged or defective instruments. Contact your local representative or dealer for repair or replacement.*

*The use of an instrument for tasks other than those for which they are intended may result in damaged / broken instruments or patient injury.*

- The operating surgeon and operating room team must be thoroughly familiar with the operating technique, as well as the range of implants and instruments to be applied. Complete information on these subjects must be readily available at the workplace.
- The operating surgeon must be especially trained in spinal surgery, biomechanical principles of the spine and the relevant operating techniques.

The operating procedure must be explained to the patient, and the patient's understanding of the following information must be documented:

- The patient is aware of the risks associated with neurosurgery, general surgery, orthopedic surgery and with general anesthesia.
- The patient has been informed about the advantages and disadvantages of the implant procedure and about possible alternative treatments.
- The patient has been informed that alternate fusion technologies (e.g., pedicle screw based systems) may be utilized, particularly if the patient's facet joints are inadequate.
- The implant can fail due to excessive load, wear and tear, or infection.
- The service life of the implant is determined by body weight and physical activity. The implant must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may be necessary if the implant fails or if there is a failure of fusion.
- The patient must have their physician carry out follow-up examinations of the implant at regular intervals.

### Intraoperative

Prior to use, please read the CoFix™ product brochure and become familiar with the surgical technique, the corresponding implants and the instruments.

- Prior to use, verify the integrity of the sterile packaging. Never use implants if the packaging is damaged.

- Prior to use, check the product expiration date. Never use implants that are past their expiration date.

The CoFix™ implant is used to provide segmental stabilization to increase the segmental stiffness and promote the fusion process. Perform microsurgical decompression or prepare central interlaminar space for CoFix™ interlaminar implant insertion, and measure interspinous space utilizing the trial portion of the coflex-F™ trial. Load the CoFix™ interlaminar implant on the insertion portion of the coflex-F™ trial and introduce into the prepared site. The wings of the CoFix™ interlaminar implant are opened with the coflex™ bending pliers and closed with the coflex™ crimping pliers. The spinous process is then prepared using the coflex-F™ punching pliers, and the screw/sleeve is inserted upon the crimped wings and fully tightened using the applicable coflex-F™ instruments.

#### Postoperative

Reiterate preoperative instructions to the patient. Ensure that the patient is aware of physical activity restrictions and possible adverse reactions.

#### Revision surgery / implant removal

The CoFix™ implant is intended for permanent implantation and is not intended for removal. However, removal may be advisable in the following situations:

- Implant breakage
- Pain due to the implant
- Infection
- Failure of fusion

#### Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use.

Technical alterations reserved.

#### For further information

Please contact Paradigm Spine or your authorized representative if further information on this product is needed.

#### Patent

For patent information see [www.surgalign.com](http://www.surgalign.com)

#### Symbols

	Manufacturer
	Date of manufacture
	Use-by date
	Lot number
	Catalog number
	Content of usable units(s)
	Sterilized using irradiation
	Do not resterilize
	Do not use if package is damaged
	Keep dry
	Do not reuse
	Consult instructions for use at this website.
	Caution
	MR conditional
	CE-marking according to directive 93/42/EEC

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