Fortilink[®]-C
with TiPlus Technology

Fortilink – C with TiPlus Technology Inserter

Reprocessing instructions

Disassembly instructions

Instructions For Use | EN

Package contents

The package contains the Fortilink - C inserter. This is either packed individually in a polyethylene bag or stored in the sterilization container as part of the instrument set.

Material

The Fortilink-C TiPlus inserter is made from the materials defined below:

- Stainless Steel per ASTM F899
- Sulfone Plastic
- Aluminum 6061-T6

Description

The Fortilink-C instruments, whilst out of the “medical device” category, are yet intended to be used specifically with the TiPlus-C implants to enable them to be applied in accordance with their intended purpose.

The surgical instruments are intended to be used for spinal instrumentation. The actual therapeutic effect is achieved by the implant. The sterilization tray is intended to hold instruments during transport, steam sterilization and storage.

The Fortilink-C Inserter is a reusable device that is provided within the dedicated Fortilink-C System steam sterilization case with other instruments for the Fortilink-C system to be sterilized by the end user. It is a device-specific accessory intended to be used with Fortilink with TiPlus Technology Interbody devices for the cervical spine. The design requirements of the Fortilink-C Inserter include:

1. Attach to interbody. Attachment is required to retain, insert, manipulate, and remove the interbody.
2. Release the interbody once positioning is complete.
3. Reattach to the interbody for device repositioning or removal if necessary.

The instruments in this system are intended for use in surgical procedures of cervical spine (C2-T1) in patients with degenerative disc diseases. These include dedicated instruments such as Fortilink C Inserters which are provided specific to each implant footprint to allow for implant insertion and removal. These device specific accessories are included in system specific steam sterilization cases with other general instruments for the Fortilink with TiPlus Technology system.

Intended use

The Fortilink-C TiPlus Inserters do not have an intended purpose according to Regulation (EU) 2017/745 Article 2 Number 1. The Fortilink-C Inserters are accessories of the Fortilink TiPlus Implants. Therefore, the intended purpose of the Fortilink TiPlus Implants covers the following aspects of the Fortilink-C TiPlus Inserter: indications, contraindications, intended patient population.

Indications

The Fortilink-C with TiPlus Technology is intended for spinal fusion procedures at one or more levels (C2-T1) in skeletally mature patients with degenerative disc disease (DDD) and instabilities at one or two contiguous levels from C2 to T1 with accompanying radicular symptoms, ruptured or herniated discs, and pseudarthrosis or failed spondylolysis. DDD is defined as neck pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies. The system is intended to be used with autograft and/or allogenic bone graft comprised of cancellous and/or

corticocancellous bone graft to facilitate fusion and is implanted via an anterior approach. The Fortilink-C Ti devices are intended to be used with supplemental fixation designed for the implanted level. This system is to be used in patients who have had six weeks of non-operative treatment.

Contraindications

The devices are contraindicated and should not be used including but not limited to, patients with the following:

- Active systemic infection or an active infection at the operative site
- A demonstrated allergy or sensitivity to any of the implant materials
- Severe osteoporosis
- Primary or metastatic tumors affecting the spine
- Conditions that may place excessive stress on bones and the implants, including but not limited to morbid obesity or other degenerative diseases
- Patients whose ability to follow postoperative restrictions, precautions and rehabilitation programs is limited
- Fractures, severe deformities or a severe instability in the area of surgery
- A medical or surgical situation that would preclude the benefit or surgery
- Pregnancy

Potential adverse effects

- Pain, discomfort, or abnormal sensations
- Dural leak
- Bone fracture
- Allergic reaction
- Cutting of skin or gloves of operating staff
- Vascular or visceral injury
- Neurological injury
- Infection
- Hemorrhage
- Death

Safety precautions

Prior to use the surgeon must become familiar with the device system and the surgical procedure. Use surgical instrumentation, accessories, and surgical technique guide provided with this device system.

- Prior to use, thoroughly read these instructions for use.
- Keep the instructions for use accessible to the staff.
- Placement and positional adjustment of implants must only be performed with system specific instruments. They must not be used with other instrumentation unless specifically recommended by Surgalign or BAAT Medical because the combination may be incompatible.
- The implants should not be rotated during insertion. Attempting to rotate these implants during insertion may lead to implant failure.
- The condition of all instruments should be checked prior to use. Damaged or worn instruments may present significant risks to safety and/or inability to function as intended and should not be used. Surgalign customer service should be contacted for instrument replacement.

Instructions For Use | EN

- If instruments are damaged or broken during use, metal fragments can be viewed by radiographic assessment. It is the surgeon's responsibility to carefully consider the risks and benefits of retrieving the fragments.
- Preoperative and operating procedures, including knowledge of surgical techniques, are important considerations in the successful utilization of the system by the surgeon. Further, the proper selection and the compliance of the patient will greatly affect the results.
- The surgeon's experience and training in surgical procedures affect the performance of this system and the use of instruments. A successful result is not always achieved in every surgical case where patient conditions may compromise the results.
- Proper patient selection and operative care are critical to the success of the surgery and avoidance of injury during surgery. Read and follow all other product information supplied by the manufacturer of the implants or the instruments.
- The surgeon should consider the levels of implantation, patient weight, patient activity level, and other patient conditions which may impact the performance of the device.
- Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without previous surgery.
- If this device is used in a patient with, or suspected with, Creutzfeldt-Jakob disease (CJD), variant CJD or other transmissible spongiform encephalopathy (TSE) and related infections, safely dispose of all the contaminated devices in accordance with local procedures and guidelines.
- The use of a surgical instrument for tasks other than those for which they are intended may result in damaged or broken instruments or patient injury.
- Follow the general guidelines and aseptic principles when handling items to be sterilized or sterile items.
- Wear appropriate protective gloves, eyewear and clothing when handling biologically contaminated instruments.
- Caution should be used while handling delicate instruments or instruments with sharp tips.
- The personnel entrusted with the processing must have the requisite knowledge and qualifications.
- Please ensure that only appropriate equipment is used for cleaning, disinfection, and sterilization of the reusable surgical devices and that performance requalification, periodic routine tests are carried out on the equipment in accordance with local procedures, guidelines, and standards.

Warnings

The same medical/ surgical conditions or complications that apply to any surgical procedure may also occur during or following implantation of this device system. The surgeon is responsible for informing the patient of the potential risks associated with treatment, including complications and adverse reactions. The surgeon may need to perform additional surgery to address any complications or adverse reactions, which may or may not be device related.

The safety, effectiveness and performance of the devices have been established for conditions in which the devices are used as intended and when used as identified in the Indications for Use. The performance of the devices has not been evaluated for use that is

contrary to the Intended Use, Indications for Use or for use that is contraindicated. Failure to use the devices as indicated could detrimentally affect the performance of their components. Breakage, slippage, misuse, unintended disassembly of multi-component instruments or mishandling of instruments, such as on sharp edges, may cause injury to the patient or operative personnel and/ or increased operating time.

Improper maintenance, handling, or poor cleaning procedures can render the instrument unsuitable for its intended purpose or even dangerous to the patient or surgical staff. These instruments should be carefully placed on trays, cleaned after each use, and stored in a dry environment.

It is important the surgeon exercises extreme caution when working near vital organs, nerves, or vessels and the forces applied while correcting the position of the instrumentation is not excessive, such that it might cause injury to the patient.

General instructions for use

This package insert is designed to assist in using surgical instruments only and is not intended to provide information on surgical technique for specific Fortilink with TiPlus Technology devices. Contact Surgalign customer service for specific product Surgical Technique Manual.

CLEANING & REPROCESSING

Fundamentals

These reusable surgical instruments are delivered unsterile. Therefore, they are to be cleaned, disinfected, and sterilized prior to each application. Effective cleaning and disinfection are indispensable requirements for an effective sterilization of the instruments.

As health care provider, additionally pay attention to the legal provisions valid for your country as well as to the hygienic requirements. This applies particularly to the different guidelines regarding the inactivation of prions (not relevant for USA), which can require the application of cleaning detergents with proven prion efficiency as well as a sterilization with more intensive parameters.

Non-Sterile Instruments:

- Instruments must be cleaned and sterilized per the below instructions prior to introduction into a sterile surgical field. or (if applicable) return of the product to the manufacturer.
- To minimize corrosion and prolong the usable life of instruments used during surgery, remove gross traces of blood and residues then thoroughly clean and dry immediately after use. Do not allow soils to dry.
- Never use steel brushes or abrasive pads, as these rupture the passive layer of the instrument surface which can lead to corrosion.
- Prior to and during use, including reprocessing, inspect instruments for:
 - Damage such as but not limited to, wear, discoloration, corrosion, cracking, fracture, or unrecognizable markings.
 - Proper function including but not limited to, sharpness, movement of hinges and couplings, joint stability, and legible markings.
 - Instruments that show signs of damage or an inability to function should not be used and should be returned to the manufacturer.

Instructions For Use | EN

Pretreatment, cleaning and disinfection

- If possible, an automated procedure (with Washer-Disinfector) should be used for cleaning and disinfection of the instruments. A manual procedure – even in case of application of an ultrasonic bath – should only be used if an automated procedure is not available; in this case, the significantly lower efficiency and reproducibility of a manual procedure has to be considered.
- It is recommended that devices should be processed as soon as reasonably possible after surgical use to facilitate the cleaning process afterwards.
- Instruments must be cleaned separately from instrument trays and cases.
- Do not soak instruments in any solution for more than two hours.
- The cleaning and disinfection detergent should be aldehyde-free (otherwise fixation of blood impurities) and should possess a fundamentally approved efficiency (for example VAH/DGHM or FDA/EPA approval/ clearance/ registration or CE marking), be suitable for the disinfection of instruments made of metallic or plastic material.
- An alkaline enzymatic detergent with low-foam formulation is recommended.
- Please consider, that a disinfectant used in the pre-treatment step serves only the personnel's safety but cannot replace the disinfection step later to be performed after cleaning.

The fundamental suitability of the instruments for an effective automated cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the Washer Disinfector G 7836 CD, Miele & Cie. GmbH & Co., Gütersloh, (thermal disinfection) and the cleaning detergent Neodisher mediclean forte (an alkaline enzymatic detergent from Dr. Weigert GmbH & Co. KG, Hamburg) considering the specified procedure.

Procedure: Pretreatment

- Multi-part or complex instruments may require disassembly for effective cleaning and disinfection. Refer additional disassembly instructions included in this IFU.
- Rinse the instruments at least 1 minute under running water (temperature < 35 °C/95 °F). If applicable actuate movable parts at least 3 times thereby. Thoroughly flush all lumens or cavities of the instruments at the beginning of the soaking time. Use thereby a syringe or a jet pistol to flush until no more soil is visible or to apply at least three times a rinsing volume of 50ml or greater.
- Soak the disassembled instruments for the recommended soaking time (Reference the instructions from the detergent manufacturer regarding concentration and soaking time) in the pre-cleaning solution so that the instruments are sufficiently covered. Pay attention not to overlap the instruments.
- Assist cleaning by careful brushing of all inner and outer surfaces with an appropriate soft brush. Brush size should be approximately the same diameter and at least the same length as of the lumen to be cleaned. Using a brush that is too big or too small in the diameter or too short may result in ineffective cleaning of the lumen/ cannulation.
- Rinse the instruments at least 1 minute under running water (temperature < 35 °C/95 °F). If applicable actuate movable parts at least 3 times during pre-rinsing. Thoroughly flush all

lumens or cavities of the instruments at the end of the soaking time. Use thereby a syringe or a jet pistol to flush until no more soil is visible or to apply at least three times a rinsing volume of 50ml or greater.

- Activate ultrasound for an additional soaking time (but not less than 5 min).
- Then, remove the instruments from the pre-cleaning solution, final rinse them intensively at least 1 minute under running water (temperature < 35 °C/95 °F). If applicable actuate movable parts at least 3 times during post-rinsing. Thoroughly flush all lumens or cavities of the instruments at the end of the soaking time. Use thereby a syringe or a jet pistol to flush until no more soil is visible or to apply at least three times a rinsing volume of 50ml or greater.
- Dry instruments with a soft, clean, absorbent non-shedding lint-free cloth. Blow clean with compressed filtered air through areas difficult to access and all lumens.

Prepare the cleaning solution/bath according to the recommendations of the detergent manufacturer regarding concentration, water quality, temperature and soaking time as well as post-rinsing and suitability for ultrasonic bath. Please use only freshly prepared solutions, only sterile or low contaminated water (max. 10 germs/ml) as well as low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water.

Procedure: Automated Cleaning and Disinfection

- Use a Washer-Disinfector (WD) with proven efficacy (validated according to EN ISO 15883 or with DGHM / FDA clearance) and regular maintenance.
- Load the disassembled, pre-treated instruments (refer additional disassembly instructions) in the WD making sure the that the lumens/cannulations can drain freely. If applicable connect the instruments to the rinsing ports of the WD. Caution shall be taken not to overlap the instruments.
- Process instruments using a validated standard WD instrument cycle in accordance with the WD manufacturer's instruction taking into account following minimal setting:

Cycle	Cycle time	Process temperature	Medium
Pre-rinsing	1 min	< 40 °C (104 °F)	Tap water
Cleaning	5 min	55 °C ± 4 (131 °F)	Alkaline enzymatic detergent
Rinsing	1 min	> 40 °C (104 °F)	Deionized water
Thermal disinfection	5 min	90 °C ± 3 (194 °F)	AD
Drying	25 min	100 °C (212 °F)	Filtered Air

- Disconnect and remove the instruments from the WD after end of the program.
- Check and pack the instruments immediately (see section "Inspection, Maintenance and Functional check") if necessary, after post-drying at a clean place).

Note:

The thermal disinfection cycle should be performed to achieve a value $A_0 \geq 3000$ (e.g. 90°C/194°F for at least 5 minutes according to ISO 15883-1).

Instructions For Use | EN

Use purified water (distilled water or de-ionized water, H₂O pure or demineralized water where the application requires this standard) for post-rinsing and thermal disinfection.

Only use of filtered air (oil-free, low contamination with microorganisms and particles) for drying.

Inspection, Maintenance and Functional check

Inspection

The cleaned and disinfected instruments are visually inspected for soil or detergent residues, damages and moisture. Worn-out, corroded, deformed or otherwise damaged instruments are not to be further reprocessed. Instruments that are still contaminated at this stage are sorted out and subjected once more to the cleaning and disinfection process. In cases of damages, the cleaned and disinfected instruments are sterilized and made available for repair.

Never use or process damaged or defective instruments. Contact your local representative or dealer for repair or replacement.

If moisture is detected, manual drying should be performed.

Maintenance

No maintenance is carried out by the processing facility.

Functional check and lubrication

Functional check should be performed where possible. Improperly functioning devices, devices with unrecognizable markings, missing or removed (buffed off) part numbers, damaged and worn devices should not be used. The disassembled instruments should be reassembled for the functional check.

Packaging

Place the cleaned, disinfected, inspected and dry devices on the specified holders within the sterilization cases provided.

Only legally marketed, and locally approved sterilization barriers (e.g. wraps, pouches or containers) should be used in compliance with the manufacturer's instructions for packaging terminally sterilized devices.

Place the sterilization case in a sterilization container, which fulfills following requirements:

- EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature resistance up to a minimum of 140 °C / 284 °F, sufficient steam permeability)
- sufficient protection of the instruments as well as of the sterilization packagings against mechanical damage
- regular maintenance of the sterilization container according to the instructions of the manufacturer

The total weight of the instrument case should not exceed 10.8kg (for metallic tray) and 6.2 kg (for plastic sterilization trays).

STERILIZATION

Independent testing has shown the following conditions to be effective:

Method	Steam
Cycle	Pre-vacuum (Wrapped)
Temperature	134°C (273°F)
Recommended Exposure Time	5-18 Minutes
Recommended Dry Time	40 Minutes

- Steam sterilization (moist heat) is the only method permitted for sterilization. Other sterilization methods are not permitted.
- The values specified here (duration/temperature) can achieve a sterility assurance level (SAL) of at least 10⁻⁶ (according to BS EN 556).
- Use a steam sterilizer with demonstrated efficacy in accordance with EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance), which is validated according to EN ISO 17665, properly maintained and calibrated. It is the responsibility of the medical facility to ensure that reprocessing is performed using the appropriate equipment and materials, and that personnel in the reprocessing facility have been adequately trained.

Do not use flash sterilization procedure, dry heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, as well as plasma sterilization.

Transportation

Surgically used instruments are safely transported to the reprocessing site in a closed container.

The use of a fractionated vacuum steam sterilization process is recommended with following validated parameters:

Cycle type	Pulses	Exposure temperature	Minimum exposure time	Minimum drying time
Dynamic air removal	3	134 °C (273 °F)	5-18 Minutes	40

Drying time for instruments processed in containers and wrapped trays can vary depending upon the type of packaging, type of instruments, type of sterilizer and total load. A minimum drying time of 30 minutes is recommended but might be extended if necessary.

Use a steam sterilizer with demonstrated efficacy in accordance with EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance), which is validated according to EN ISO 17665, properly maintained and calibrated. It is the responsibility of the medical facility to ensure that reprocessing is performed using the appropriate equipment and materials, and that personnel in the reprocessing facility have been adequately trained.

Do not use flash sterilization procedure, dry heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, as well as plasma sterilization.

Material resistance

The following substances must not be ingredients of the cleaning detergent:

- organic, mineral, and oxidizing acids (minimum admitted pH-value 5.5)
- strong lyes (maximum admitted pH-value 11, weak alkaline cleaner with prion effectiveness recommended)
- organic solvents (for example: acetone, ether, alcohol, benzine)
- oxidizing agents (for example: peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic, halogenated hydrocarbons

It is not recommended to add corrosion inhibitors, neutralizing agents, and/or rinse aids to the detergents; they may potentially cause critical remnants on the instruments.

Instructions For Use | EN

Do not expose any instruments or sterilization trays to temperatures equal or greater than 140°C/284°F.

Limitations on reprocessing

Repeated processing cycles in compliance with these instructions for use have minimal effects on device life cycle and function. Instruments do not have an indefinite life cycle. Their end of life is determined by wear and damage due to surgical use and handling. Evidence of damage and wear on a device may include but is not limited to corrosion (i.e. rust, pitting), discoloration, excessive scratches, flaking, wear and cracks. Improperly functioning devices, devices with unrecognizable markings, missing or removed (buffed off) part numbers, damaged and excessively worn devices should not be further used.

A maximum of 100 reprocessing cycles is recommended.

Storage

After sterilization, the reusable instruments should be stored in dry and dust-free hospital environmental conditions.

Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use. Technical alterations reserved.

Requirements regarding the information supplied with the device

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

For further information

Please contact Paradigm Spine or your authorized representative if further information on this product is needed.

Symbols

	Manufacturer
	Lot number
	Catalog number
	Country of manufacture
	Non-sterile
	Do not use if package is damaged
	Keep dry
	Consult instructions for use at this website.
	Caution
	Content of usable units(s)

	Medical device
	CE-marking according to Regulation (EU) 2017/745 Annex V
	Federal law restricts this device to sale by or on the order of a physician.

Disassembly Instructions


Do not use a tool with sharp ends.

Fortlink-C TiPlus Inserter 65-C-INS-TIPLUS

Step 1: Compress inserter fork distal tip and rotate proximal handle.

Step 2: Inserter fork will translate from the inserter housing until disassembly is complete.

Fortilink®-C
with TiPlus Technology

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