
Instructions for Use

Reusable surgical instruments and sterilization trays

Reprocessing instructions

Disassembly instructions

Reusable surgical instruments and sterilization trays

Package contents

The package contains one equipped sterilization tray.

Description

The sterilization tray is a reusable tray that holds the surgical instruments before, during and after steam sterilization.

The Paradigm Spine reusable surgical instruments shipped with this sterilization tray are designed for use in orthopedic surgical procedures for implantation of devices manufactured by Paradigm Spine. Prior to use of these instruments, please refer to the implant instructions for use specific to the product line and the surgical procedure to be followed.

Material

The sterilization tray is made from stainless steel and Polyphenylsulfone (PPSU) and may include Polypropylene (PP) brackets and/or Silicone.

The surgical instruments are mainly made from medical grade acetal copolymer (POM), stainless steel and may include silicone, aluminium alloy, titanium alloy, Polyamide 6 and/or Polytetrafluoroethylene (PTFE).

Intended use

The sterilization tray is intended to hold instruments during transport, steam sterilization and storage.

The surgical instruments are intended to be used for spinal instrumentation. The actual therapeutic effect is achieved by the implant.

Indications

Please see the IFU of the specific implant that may be used with these surgical instruments. These surgical instruments are designed to manipulate tissue or for use with other devices in orthopedic and spine surgery.

Contraindications

Please see the IFU of the specific implant that may be used with these surgical instruments. With any surgery, the following contraindications exist:

- Infection in or around the operative site
- Allergy or sensitivity to instrument materials
- Use of incompatible materials from other systems
- Any case not described in the indication

Adverse reactions

Please see the IFU of the specific implant that may be used with these surgical instruments.

Safety precautions

- Prior to use, thoroughly read these instructions for use.
- Keep the instructions for use accessible to the staff.
- The methods of use of these instruments are to be determined by the user's experience and training
- Never use or process damaged or defective devices. Contact your local sales representative or Paradigm Spine for repair or replacement.
- If this device is used in a patient with, or suspected with, Creutzfeldt-Jakob disease (CJD), variant CJD or other transmissible spongiform encephalopathy (TSE) and related infections, safely dispose and/or incinerate all the contaminated devices in accordance with local procedures and guidelines.
- The use of a surgical instrument for tasks other than those for which they are intended may result in damaged or broken instruments or patient injury.

- Follow the general guidelines and aseptic principles when handling items to be sterilized or sterile items.
- Wear appropriate protective gloves, eyewear and clothing when handling biologically contaminated instruments.
- Caution should be used while handling delicate instruments or instruments with sharp tips.
- Use only soft brushes, no metal brushes.
- Do not process the surgical instruments at a temperature equal or greater than 140°C/284°F.
- The personnel entrusted with the processing must have the requisite knowledge and qualifications.
- Please ensure that only appropriate equipment is used for cleaning, disinfection, and sterilization of the reusable surgical devices and that performance requalification, periodic routine tests are carried out on the equipment in accordance with local procedures, guidelines, and standards.

Processing and reprocessing of reusable surgical devices

Fundamentals

Paradigm Spine's reusable surgical instruments are delivered unsterile. Therefore, they are to be cleaned, disinfected, and sterilized prior to each application. Effective cleaning and disinfection are indispensable requirements for an effective sterilization of the instruments.

As health care provider, additionally pay attention to the legal provisions valid for your country as well as to the hygienic requirements. This applies particularly to the different guidelines regarding the inactivation of prions (not relevant for USA), which can require the application of cleaning detergents with proven prion efficiency as well as a sterilization with more intensive parameters.

Instrument background

Surgical instruments are designed to perform specific functions such as cutting, drilling, sawing, scratching, scraping, clamping, retracting, clipping or similar procedures. Surgical instruments are also used to facilitate the insertion of surgical implants. Proper cleaning, disinfection, sterilization, handling and routine maintenance help to ensure that surgical instruments perform as intended.

Paradigm Spine's surgical instruments are supplied non-sterile and must be thoroughly cleaned, disinfected and sterilized prior to use according to the instructions in this document.

The equipped sterilization tray is supplied individually packed in protective packaging that is labeled according to its contents. Before processing, inspect the instruments and sterilization tray carefully for damage, wear and/or non-functioning parts. Carefully inspect the inaccessible areas, joints and all moveable parts.

Care and handling at the point of use

As soon as possible after use, dismantle where applicable the instruments, rinse and flush lumened devices. Only disposable wipes and soft brushes (e.g. brushes with nylon bristles) may be used to remove coarse impurities.

Transportation

Surgically used instruments are safely transported to the reprocessing site in a closed container.

Pretreatment, cleaning and disinfection

- An automated procedure (with Washer-Disinfector) should be used for cleaning and disinfection of the instruments.

- It is recommended that devices should be processed as soon as reasonably possible after surgical use to facilitate the cleaning process afterwards.
- Instruments must be cleaned separately from instrument trays and cases.
- Do not soak instruments in any solution for more than two hours.
- The cleaning and disinfection detergent should be aldehyde-free (otherwise fixation of blood impurities) and should possess a fundamentally approved efficiency (for example VAH/DGHM or FDA/EPA approval/ clearance/ registration or CE marking), be suitable for the disinfection of instruments made of metallic or plastic material.
- An alkaline enzymatic detergent with low-foam formulation is recommended.
- Please consider, that a disinfectant used in the pre-treatment step serves only the personnel's safety but cannot replace the disinfection step later to be performed after cleaning.

The fundamental suitability of the instruments for an effective automated cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the Washer Disinfector G 7836 CD, Miele & Cie. GmbH & Co., Gütersloh, (thermal disinfection) and the cleaning detergent Neodisher mediclean forte (an alkaline enzymatic detergent from Dr. Weigert GmbH & Co. KG, Hamburg) considering the specified procedure. Sonorex super RK 514 H, Bandelin electronic was used as ultrasonic device with a frequency of 35 kHz.

Procedure: Pretreatment

- Multi-part or complex instruments may require disassembly for effective cleaning and disinfection. Refer additional disassembly instructions included in this IFU.
- Rinse the instruments at least 1 minute under running water (temperature < 35 °C/95 °F). If applicable actuate movable parts at least 3 times thereby. Thoroughly flush all lumens or cavities of the instruments at the beginning of the soaking time. Use thereby a syringe or a jet pistol to flush until no more soil is visible or to apply at least three times a rinsing volume of 50ml or greater.
- Soak the disassembled instruments for the recommended soaking time (Reference the instructions from the detergent manufacturer regarding concentration and soaking time) in the pre-cleaning solution so that the instruments are sufficiently covered. Pay attention not to overlap the instruments.
- Assist cleaning by careful brushing of all inner and outer surfaces with an appropriate soft brush. Brush size should be approximately the same diameter and at least the same length as of the lumen to be cleaned. Using a brush that is too big or too small in the diameter or too short may result in ineffective cleaning of the lumen/ cannulation.
- Rinse the instruments at least 1 minute under running water (temperature < 35 °C/95 °F). If applicable actuate movable parts at least 3 times during pre-rinsing. Thoroughly flush all lumens or cavities of the instruments at the end of the soaking time. Use thereby a syringe or a jet pistol to flush until no more soil is visible or to apply at least three times a rinsing volume of 50ml or greater.
- Activate ultrasound for an additional soaking time (but not less than 5 min).
- Then, remove the instruments from the pre-cleaning solution, final rinse them intensively at least 1 minute under running water (temperature < 35 °C/95 °F). If applicable actuate movable parts at least 3 times during post-rinsing. Thoroughly flush all lumens or

cavities of the instruments at the end of the soaking time. Use thereby a syringe or a jet pistol to flush until no more soil is visible or to apply at least three times a rinsing volume of 50ml or greater.

- Dry instruments with a soft, clean, absorbent non-shedding lint-free cloth. Blow clean with compressed filtered air through areas difficult to access and all lumens.

Prepare the cleaning solution/bath according to the recommendations of the detergent manufacturer regarding concentration, water quality, temperature and soaking time as well as post-rinsing and suitability for ultrasonic bath. Please use only freshly prepared solutions, only sterile or low contaminated water (max. 10 germs/ml) as well as low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water.

Procedure: Automated Cleaning and Disinfection

- Use a Washer-Disinfector (WD) with proven efficacy (validated according to EN ISO 15883 or with DGHM / FDA clearance) and regular maintenance.
- Load the disassembled, pre-treated instruments (refer additional disassembly instructions) in the WD making sure the that the lumens/cannulations can drain freely by choosing a suitable orientation and fixating them. If applicable connect the instruments to the rinsing ports of the WD. Caution shall be taken not to overlap the instruments.
- Process instruments using a validated standard WD instrument cycle in accordance with the WD manufacturer's instruction taking into account following minimal setting:

Cycle	Cycle time	Process temperature	Medium
Pre-rinsing	1 min	< 40 °C (104 °F)	Tap water
Cleaning	5 min	55 °C ± 4 (131 °F)	Alkaline enzymatic detergent
Rinsing	1 min	> 40 °C (104 °F)	Deionized water
Thermal disinfection	5 min	90 °C ± 3 (194 °F)	AD
Drying	25 min	100 °C (212 °F)	Filtered Air

- Disconnect and remove the instruments from the WD after end of the program.
- Check and pack the instruments immediately (see section "Inspection, Maintenance and Functional check") if necessary, after post-drying at a clean place).

Note:

The thermal disinfection cycle should be performed to achieve a value $A_0 \geq 3000$ (e.g. 90°C/194°F for at least 5 minutes according to ISO 15883-1).

Use purified water (distilled water or de-ionized water, H₂O pure or demineralized water where the application requires this standard) for post-rinsing and thermal disinfection.

Only use of filtered air (oil-free, low contamination with microorganisms and particles) for drying.

Inspection, Maintenance and Functional check

Inspection

The cleaned and disinfected instruments are visually inspected for soil or detergent residues, damages and moisture. Worn-out, corroded, deformed or otherwise damaged instruments are not to be further reprocessed. Instruments that are still contaminated at this stage are sorted out and subjected once more to the cleaning and disinfection process. In cases of damages, the cleaned and disinfected instruments are sterilized and made available for repair.

Never use or process damaged or defective instruments. Contact your local representative or dealer for repair or replacement.

If moisture is detected, manual drying should be performed.

Maintenance

No maintenance is carried out by the processing facility.

Functional check and lubrication

Functional check should be performed where possible. Improperly functioning devices, devices with unrecognizable markings, missing or removed (buffed off) part numbers, damaged and worn devices should not be used. The disassembled instruments should be reassembled for the functional check.

Use only instrument oils (white oil) admitted to steam sterilization considering the maximum possible sterilization temperature and with approved biocompatibility. Apply only a small amount and only to the joints.

Packaging

Place the cleaned, disinfected, inspected and dry devices on the specified holders within the sterilization cases provided. In case the single devices, which are delivered optionally does not fit in the regular sterilization case, please use another suitable medical grade sterilization pouch.

Only legally marketed, and locally approved sterilization barriers (e.g. wraps, pouches or containers) should be used in compliance with the manufacturer's instructions for packaging terminally sterilized devices. Place the sterilization case in a sterilization container, which fulfills following requirements:

- EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature resistance up to a minimum of 140 °C / 284 °F, sufficient steam permeability)
- sufficient protection of the instruments as well as of the sterilization packagings against mechanical damage
- regular maintenance of the sterilization container according to the instructions of the manufacturer

The total weight of the instrument case should not exceed 10.8kg (for metallic tray) and 6.2 kg (for plastic sterilization trays).

Sterilization

The use of a fractionated vacuum steam sterilization process is recommended with following validated parameters:

• In Germany

Cycle type	Pulses	Exposure temperature	Minimum exposure time	Minimum drying time
Dynamic air removal	3	134 °C (273 °F)	5 min	30

• In USA

Cycle type	Pulses	Exposure temperature	Minimum exposure time	Minimum drying time
Dynamic air removal	3	132 °C (270 °F)	4 min	30

• In other countries

Cycle type	Pulses	Exposure temperature	Minimum exposure time	Minimum drying time
Dynamic air removal	3	132 °C- 134 °C (270 °F-273 °F)	4 min	30
		134 °C (273 °F)	18 min ¹	

The fundamental suitability of the instruments for an effective steam sterilization was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the steam sterilizer Selectomat PL 669 (MMM Sterilisatoren AG) and the fractionated vacuum/dynamic air removal procedure and the instrument oil Sterilit i (JG600), Aesculap.

Sterilization containers with paper filter and double pouches according to EN ISO 11607 (for individually packed instruments) were used.

Drying time for instruments processed in containers and wrapped trays can vary depending upon the type of packaging, type of instruments, type of sterilizer and total load. A minimum drying time of 30 minutes is recommended but might be extended if necessary.

Use a steam sterilizer with demonstrated efficacy in accordance with EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance), which is validated according to EN ISO 17665, properly maintained and calibrated. It is the responsibility of the medical facility to ensure that reprocessing is performed using the appropriate equipment and materials, and that personnel in the reprocessing facility have been adequately trained.

Do not use flash sterilization procedure, dry heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, as well as plasma sterilization.

Material resistance

The following substances must not be ingredients of the cleaning detergent:

- organic, mineral, and oxidizing acids (minimum admitted pH-value 5.5)
- strong lyes (maximum admitted pH-value 11, weak alkaline cleaner with prion effectiveness recommended)
- organic solvents (for example: acetone, ether, alcohol, benzene)
- oxidizing agents (for example: peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic, halogenated hydrocarbons

¹ For the inactivation of prions.

It is not recommended to add corrosion inhibitors, neutralizing agents, and/or rinse aids to the detergents; they may potentially cause critical remnants on the instruments.

Do not expose any instruments or sterilization trays to temperatures equal or greater than 140°C/284°F.

Limitations on reprocessing

Repeated processing cycles in compliance with these instructions for use have minimal effects on device life cycle and function. Instruments do not have an indefinite life cycle. Their end of life is determined by wear and damage due to surgical use and handling.

Evidence of damage and wear on a device may include but is not limited to corrosion (i.e. rust, pitting), discoloration, excessive scratches, flaking, wear and cracks. Improperly functioning devices, devices with unrecognizable markings, missing or removed (buffed off) part numbers, damaged and excessively worn devices should not be further used.

A maximum of 100 reprocessing cycles is recommended.

Storage

After sterilization, the reusable instruments should be stored in dry and dust-free hospital environmental conditions.

Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use. Technical alterations reserved.

For further information

Please contact Paradigm Spine or your authorized representative if further information on this product is needed.

Symbols

	Manufacturer
	Lot number
	Catalog number
	Content of usable units(s)
	Do not use if package is damaged
	Keep dry
	Non-sterile
	Consult instructions for use at this website.
	Caution
	CE-marking according to directive 93/42/EEC
	CE-marking according to directive 93/42/EEC
	Medical Device

Disassembly instructions for specific instruments


If necessary, please use a surgical tool to press the buttons necessary to enable disassembly; for instance, a surgical probe.

Do not use a tool with sharp ends.


DCI [™] reusable devices

REF CBT20100 or CBT20110 Insertion instrument (left)

(A) dismantle by screwing of the safety stop shaft in the middle towards the front

REF HAT90900 Stiffness gauge


DSS [®] reusable devices

REF DAT40105 Pedicle screw inserter

(A) Press the notch

(B) Turn the handle till it slides down the shaft


(C) Lift the window

REF HAT5026K (HAT50260+HAT50005) Cannulated pedicle awl w/ K-wire. See DAT5056K.

REF DAT5026K (DAT50260+DAT50005) Cannulated pedicle awl w/ K-wire

1. Unscrew the K-wire with cap
2. Clean the inner and outer threads
3. Sterilize K-wire separately from awl


HPS [™] reusable devices

REF HAT90500 Persuader

(A) Set the persuader in neutral position by pushing the trigger



(B) Press the notch

coflex-F® reusable devices



REF RAT20215 Screw inserter

1,2,3 Press the buttons with the angled hook, lift the trigger and let the cover slide
4,5, Press the buttons with the angled hook and let the cover slide.

REF RAT20222 LIS screw inserter

1,2,3 Press the buttons with the angled hook, lift the trigger and let the cover slide

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