Instructions for Use

DSS™ Stabilization System

Paradigm Spine GmbH | Eisenbahnstrasse 84 | 78573 Wurmlingen | Germany
Phone +49 (0) 7461-9635990 | Telefax: +49 (0) 7461-963599-20
info@surgalign.com | www.surgalign.com

LBL 0656 Rev. F 2023-5

PARADIGM SPINE
the movement in spine care

DSS™ Stabilization System

Package contents

The package contains one or more of the various components of the DSS™ Stabilization System:

- DSS™ Pedicle screws
- DSS™ Dynamic coupler
- DSS™ Fusion coupler
- DSS™ Closure set

The content of each package is evident from the respective product label.

Description

The DSS™ Stabilization System is a pedicle screw based system for mono- or multisegmental posterior dynamic or rigid stabilization of the lumbosacral spine.

Material

The DSS™ Stabilization System implant components are made from titanium 6-aluminum 4-vanadium alloy (ISO 5832-3).

The used material contains no nickel.
Not made with natural rubber latex.

Intended Purpose

The DSS™ Stabilization System is a pedicle screw based implant system for posterior dynamic stabilization and/or for posterior rigid stabilization of the spine at one to three contiguous spinal segments from L1 to S1 in skeletally mature patients. It is intended for long-term implantation and for single-use only. It is delivered sterile and is not intended for re-processing by the user or third parties. It shall only be used by surgeons who are trained and familiar with the implant components, instruments, and surgical technique.

Indications

Posterior Dynamic Stabilization

The DSS™ as a system for posterior dynamic stabilization is intended for posterior dynamic stabilization after (micro-) surgical decompression in the treatment of degenerative spinal stenosis including degenerative spondylolisthesis (pseudospondylolisthesis) up to Meyerding grade 1, which were previously unsuccessfully conservatively treated.

Posterior Rigid Stabilization

The DSS™ as a posterior rigid stabilization system is intended for posterior rigid stabilization as adjunct to fusion for the following indications:

- degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies)
- spondylolisthesis
- trauma (i.e., fracture or dislocation)

Contraindications

Generally, do not use the DSS™ Stabilization System for the following conditions:

- All medical and surgical conditions that exclude the advantages of spine surgery
- Damage to the bone structure that prevents a stable implantation
- Acute or chronic systemic, spinal or localized infections
- Impaired wound healing
- Active systemic and metabolic diseases
- Obesity (BMI ≥ 30)
- Fever
- Pregnancy
- Medication and substance abuse, alcoholism
- Lack of patient cooperation

- Mental illness
- Sensitivity to foreign bodies such as implant materials
- Serious osteopenia
- Osteoporosis
- Neuromuscular disorders or illnesses
- Bone tumors in the proximity of the implant fixation
- Other cases that are not listed among the "Indications"

Do not use the DSS™ Dynamic Stabilization System in cases of:

- Degenerative spondylolisthesis greater than grade 1
- Isthmic spondylolisthesis
- Dysplastic spondylolisthesis
- Severe instabilities
- Degenerative scoliosis greater than 15 degrees
- Fracture
- Tumor
- Infection

Adverse reactions

Adverse reactions may include:

- Clinical failure (i.e. pain or injury) due to bending, loosening, wear and tear, fracture of implant, loss of fixation, dislocation and/or migration of the implant
- Pain and/or abnormal sensations due to the presence of the implant
- Primary or secondary infections
- Allergic reactions to implant material
- Neurological injury
- Vertebrae fracture
- Injury to vessels, nerves, the spinal cord, and organs
- Hematoma and/or impaired wound healing; hemorrhage
- Venous thrombosis, lung embolism, and cardiac arrest
- Lack of fixation
- No or delayed fusion (only applies to fusion components)
- Change to the spinal column curvature
- Impairment of the gastrointestinal, urological and/or reproductive system
- Pain or malaise
- Bursitis
- Decrease of bone density due to avoiding load
- Osteoporosis or fracture above or below the affected area of the spine
- Persistence of the symptoms to be treated with the implantation
- Death

Safety precautions

- Prior to use, thoroughly read these instructions for use and become familiar with the surgical technique.

- Only stabilize a segment with either dynamic couplers or fusion couplers. **Do not use both a dynamic coupler and a fusion coupler at the same level.**
- Keep the instructions for use accessible to all staff.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques. Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrect operating techniques, the limitations of treatment methods or inadequate asepsis.

- Under no circumstances may modular implant components from different suppliers be combined.
- Each patient's record shall document the implant used (name, article number, lot number).
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed about post-surgical behavioral requirements.
- Damage to the weight-bearing structures can give rise to loosening, dislocation and migration, as well as other complications.
To ensure the earliest possible detection of such catalysts of implant malfunction, the implant must be checked periodically post-operatively using appropriate techniques.

- Never reuse an implant. Although the implant may appear undamaged, previous stresses may have created non-visible damage that could result in implant failure.

- Never use implants if the packaging is damaged. An implant with damaged packaging might be damaged itself and thus may not be used.

- Never use implants that are past their expiration date.

MRI compatibility

The DSS™ stabilization system has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the DSS™ implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Storage, inspection and sterilization

Storage

The DSS™ implants are packed in protective packaging that is labeled according to its contents.

The implant is sterilized with gamma sterilization (25 kGy minimum).

- Always store the implant in the original protective packaging.
- Do not remove the implant from the packaging until immediately before use.
- Store the implant in a dry and dust-free place (standard hospital environment).

Disinfection / cleaning

The DSS™ implants are not designed to be disinfected or cleaned by the user.

Restoration

The DSS™ implants are not designed to be re-sterilized by the user.

Procedure

The DSS™ implants must be implanted only with the applicable DSS™ instruments. The DSS™ instruments are available from the manufacturer at any time.

Preoperative

The operating surgeon draws up an operation plan specifying and documenting the following:

- Implant component(s) and their dimensions.
- Proper positioning of the implant components.
- Determination of intra-operative orientation points.

The following conditions must be fulfilled prior to application:

- All required implant component(s) are readily available.

- Highly aseptic operating conditions are present.
- All requisite implantation instruments must be available and in working order.

Never use or process damaged or defective instruments. Contact your local representative or dealer for repair or replacement.

The use of an instrument for tasks other than those for which they are intended may result in damaged / broken instruments or patient injury.

- The operating surgeon must be especially trained in spinal surgery, biomechanical principles of the spine and the relevant operating techniques.
- The operating surgeon and operating room team must be thoroughly familiar with the operating technique, as well as the range of implants and instruments to be applied. Complete information on these subjects must be readily available at the workplace.
- The operating procedure must be explained to the patient, and the patient's understanding of the following information must be documented:
 - The patient is aware of the risks associated with neurosurgery, general surgery, orthopedic surgery and with general anesthesia.
 - The patient has been informed about the advantages and disadvantages of the implant procedure and about possible alternative treatments.
 - The implant can fail due to excessive load, wear and tear, or infection.
 - The service life of the implant is determined by body weight and physical activity. The implant must not be subjected to overload through extreme strain, or through work-related or athletic activities.
 - Corrective surgery may be necessary if the implant fails.
- The patient must have their physician carry out follow-up examinations of the implant at regular intervals.

Intraoperative

Prior to use, please read the DSS™ product brochure and become familiar with the surgical technique, the corresponding implants and the instruments.

- Prior to use, verify the integrity of the sterile packaging. Never use implants if the packaging is damaged.

- Prior to use, check the product expiration date. Never use implants that are past their expiration date.

Prepare the surgical site for implant introduction. Numerous instruments (e.g., trocar, awl, pedicle finder, reamer, tap and pedicle probe) are provided to prepare the pedicles prior to screw insertion.

Introduce the pedicle screws into the pedicles of the affected motion segment using the pedicle screw inserter (note: the pedicle screws are cannulated and may be used with 1.6 mm K-wires). Attach the pedicle screw extension rods to the pedicle screws (after removing K-wires, if applicable). Measure the distance between the pedicle screws using the caliper set, and select the applicable size coupler based on the measured distance. Adjust the dynamic coupler to the measured distance using the dynamic coupler length adjuster. Do not apply excessive force when adjusting the coupler length. The couplers have built-in stops that limit the minimum and maximum adjustment range.

For single level dynamic stabilization (non-fusion), apply the dynamic coupler to the pedicle screws in the orientation as indicated on the implant (the cranial direction is marked with UP and an arrow on the dynamic coupler). The coupler pusher may be used to firmly seat the dynamic coupler onto the pedicle screws.

Add a spacer using the spacer inserter, then a washer and then a nut using the nut driver to each pedicle screw (note: alternative configurations may be used by the surgeon in order to optimize the fit to the anatomy such as screw, washer, spacer, dynamic coupler, nut). Fix the dynamic coupler to the pedicle screws by tightening the nuts utilizing a counter torque instrument to avoid bending of the coupler and movement of the pedicle screws during tightening. If desired, remove the extension rods with the extension rod remover prior to applying the counter torque instrument. The proper tightening torque is reached when the upper half of the nut breaks off. The slots of the coupler shall remain parallel during tightening (note: if the slots of the coupler are not parallel, correct the coupler position using a counter torque instrument). Remove all instruments.

For single level rigid stabilization (fusion), apply the spacer inserter, then a washer and then a nut using the nut driver to each pedicle screw (note: alternative configurations may be used by the surgeon in order to optimize the fit to the anatomy such as pedicle screw, washer, spacer, fusion coupler, nut). Use a distraction and/or compression instrument to achieve sagittal alignment. Once aligned, use the hex driver to "lock" the fusion coupler. Fix the fusion coupler to the pedicle screws by tightening the nuts. A counter torque instrument shall be utilized to avoid movement of the pedicle screws. If desired, the tightening of the nuts may be done prior to "locking" the fusion coupler. Again, the proper tightening torque is reached when the upper half of the nut breaks off.

For multi-level cases, only the pedicle screws at both ends of the instrumentation receive a washer. The screws in between receive a first coupler, a spacer, a second coupler and the nut. Any additional washers contained in the blister package shall be disposed. Other than this, the proper technique is applicable as described above for rigid or dynamic stabilization.

Postoperative

- Reiterate preoperative instructions to the patient.
- Ensure that the patient is aware of physical activity restrictions and possible adverse reactions.

Revision surgery / implant removal

The DSS dynamic stabilization system provides dynamic stability to a degenerated and destabilized segment. The lifetime of the implant is 48 months. By reaching the end of the lifetime the DSS dynamic stabilization system shall be explanted and may be replaced by the DSS rigid stabilization system, as deemed necessary by the physician.

The DSS rigid stabilization system provides initial primary stability to the treated segment until fusion occurs which creates long-term secondary stability. The time to fusion reported in the literature is ≤ 24 months which is considered the expected functional lifetime of the device. The DSS rigid stabilization system shall be explanted when a solid fusion is evident.

However, removal of the implant can be necessary in the following situations:

- Implant breakage
- Pain due to the implant
- Infection
- Pseudarthrosis
- Allergic reactions

If a revision is performed where dynamic couplers have been utilized and the surgeon elects to fuse this segment, the dynamic couplers may be replaced with fusion couplers. If it is necessary to replace the dynamic coupler and/or fusion coupler, it is advisable to

replace them on both sides of the spinal segment. Implants which appear to be intact may have invisible damage. The surgeon shall determine whether new pedicle screws are used, or the coupler is placed over the existing screws. The surgeon must however use new nuts, washers (as applicable) and spacers.

To remove a DSS construct, carry out the following steps as described. Pay attention to the loosened components during the procedure.

Remove the nuts:

Use the nut driver in combination with the ratchet to remove the nuts counterclockwise. For initial release utilizing the counter torque is strongly recommended.

Removal of couplers, spacers and washers:

After removing all nuts, the couplers, washers and spacers can be removed. Care has to be taken that all these components are removed.

Screw removal:

The T-Handle is connected to the pedicle screw inserter and then inserted deeply over the top of the screws. After ensuring a good connection between screw and pedicle screw inserter the screws can be removed by turning the t-handle on the pedicle screw inserter counterclockwise.

Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with these instructions for use. Subject to technical modifications.

Disposal

Any implant that has been contaminated by blood, tissue, and/or bodily fluids can potentially cause infection and injury and should therefore never be used again. Contaminated implants should be handled according to hospital protocol.

Devices must be disposed of as a healthcare medical device in compliance with hospital procedures and national and legal requirements.

For further information

Please contact Paradigm Spine or your authorized representative if further information on this product is needed


IMPORTANT NOTE

These instructions for use are complete at the time of publishing.

Patent

For patent information see www.surgalign.com

Symbols

	Manufacturer
	Date & Country of manufacture
	Use-by date
	Batch code
	Catalogue number
	Content of usable units(s)
	Sterilized using irradiation
	Do not resterilize
	Do not use if package is damaged
	Keep dry
	Do not reuse
	Consult instructions for use at this website.
	Medical device
	Caution
	CE-marking according to directive 93/42/EEC
	Double sterile barrier system
	Unique device identifier

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