# Instructions for Use

HPS™

Hybrid Performance System

Posterior Dynamic Stabilization

Paradigm Spine GmbH | Eisenbahnstrasse 84 | 78573 Wurmlingen | Germany  
Phone +49 (0) 7461-9635990 | Telefax: +49 (0) 7461-963599-20  
[info@surgalign.com](mailto:info@surgalign.com) | [www.surgalign.com](http://www.surgalign.com)

LBL 0702 Rev. M 2023-05

**PARADIGM SPINE**  
*the movement in spine care*

## HPS™ Hybrid Performance System

### Package contents

The package contains one or more of the various components of the HPS™ Hybrid Performance System or HPS™ 2.0 Hybrid Performance System:

- HPS™ Closure set
- HPS™ 2.0 Closure set
- HPS™ Polyaxial pedicle screw
- HPS™ Polyaxial pedicle screw, Ti-coated
- HPS™ Polyaxial pedicle screw, fenestrated
- HPS™ Monoaxial pedicle screw
- HPS™ 2.0 Polyaxial pedicle screw
- HPS™ Rod coupler, S=25 N/mm
- HPS™ Rod coupler, S=50 N/mm
- HPS™ 2.0 Rod coupler, S=25 N/mm
- HPS™ 2.0 Rod coupler, S=50 N/mm
- HPS™ Rod, pre-bent
- HPS™ Rod, straight

The content of each package is evident from the respective product label.

### Description

The modular HPS™ Hybrid Performance System is used for single or multi-segmental posterior stabilization of the thoracolumbar spine. The system is based on fixation with mono- and/or polyaxial pedicle screws. The system allows both dynamic and rigid stabilization. The HPS™ Pedicle screw, fenestrated is also used for single or multi-segmental posterior rigid stabilization of the thoracolumbar spine in combination with bone cement.

### Material

The HPS™ Hybrid Performance System is made from Titanium 6-aluminum 4-vanadium alloy (ISO 5832-3).

The used material contains no nickel.

Not made with natural rubber latex.

### Intended use

The HPS™ Hybrid Performance System is a pedicle screw-based implant system for posterior dynamic stabilization and/or for posterior rigid stabilization of the thoracolumbar spine in skeletally mature patients. The HPS is intended for long-term implantation and for single-use only. It is delivered sterile and is not intended for re-processing by the user or third parties. The HPS™ Hybrid Performance System shall only be used by surgeons who are trained and familiar with the implant components, instruments, and surgical technique.

### Indications

#### Posterior Dynamic Stabilization

The HPS™ Hybrid Performance System as a system for posterior dynamic stabilization is intended for posterior dynamic stabilization at one to three contiguous spinal segments from L1 to S1 after (micro-) surgical decompression in the treatment of degenerative spinal stenosis including degenerative spondylolisthesis (pseudospondylolisthesis) up to Meyerding grade 1, which were previously unsuccessfully conservatively treated.

#### Posterior Rigid Stabilization

The HPS™ Hybrid Performance System as a posterior rigid stabilization system is intended for posterior rigid stabilization of the thoracolumbar spine as adjunct to fusion for the following indications:

- degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies);

- spondylolisthesis;
- trauma (i.e., fracture or dislocation);
- spinal stenosis;
- curvatures (i.e., kyphosis, and/or lordosis);
- tumor,
- pseudoarthrosis;
- and failed previous fusion

#### HPS™ Pedicle screw, fenestrated

The HPS™ Pedicle screw, fenestrated, is indicated for the use with bone cement in patients with diminished bone quality (e.g. osteoporosis, osteopenia, metastatic disease) and shall only be used with components for posterior rigid stabilization of the HPS™ Hybrid Performance System for mono- or multisegmental stabilization of the thoracic and lumbar spine for supporting fusion.

### Contraindications

#### Posterior Dynamic Stabilization

Do not use the HPS™ Rod coupler in cases of:

- Degenerative spondylolisthesis greater than grade 1 or isthmic spondylolisthesis
- Degenerative scoliosis greater than 15 degrees
- Dysplastic spondylolisthesis
- Fracture
- Tumor
- Infection.
- Severe instabilities

### General contraindications

Generally do not use the HPS™ Hybrid Performance System for the following conditions:

- All medical and surgical conditions that exclude the advantages of spine surgery
- Damage to the bone structure that prevents a stable implantation
- Acute or chronic systemic, spinal or localized infections
- Impaired wound healing
- Active systemic and metabolic diseases
- Obesity (BMI ≥30)
- Fever
- Pregnancy
- Medication and substance abuse, alcoholism
- Lack of patient cooperation
- Mental illness
- Sensitivity to foreign bodies such as implant materials
- Serious osteopenia
- Osteoporosis
- Neuromuscular disorders or illnesses
- Bone tumors in the proximity of the implant fixation
- Fracture
- Other cases that are not listed among the "Indications"

### Contraindications with use of bone cement

- Poor visualizing at Fluoroscopy
- Thrombophilia
- Patients with severe heart and/or lung failure
- Patients with known sensitiveness to elements of the bone cement

Refer to the IFU of the bone cement supplier for further contraindications

## Adverse reactions

Adverse reactions may include:

- Clinical failure (i.e. pain or injury) due to bending, loosening, wear and tear, fracture of implant, loss of fixation, dislocation and/or migration of the implant
- Pain and/or abnormal sensations due to the presence of the implant
- Primary or secondary infections
- Allergic reactions to implant material
- Neurological injury
- Vertebrae fracture
- Injury to vessels, nerves, the spinal cord, and organs
- Hematoma and/or impaired wound healing; hemorrhage
- Venous thrombosis, lung embolism, and cardiac arrest
- Lack of fixation
- No or delayed fusion (only applies to fusion components)
- Change to the spinal column curvature
- Impairment of the gastrointestinal, urological and/or reproductive system
- Pain or malaise
- Bursitis
- Decrease of bone density due to avoiding load
- Osteoporosis or fracture above or below the affected area of the spine
- Persistence of the symptoms to be treated with the implantation
- death

## Adverse reaction with use of bone cement

- Hypersensitivity reactions by susceptible patients which may cause to an anaphylactic shock
- Tissue damage, nerve damage or circulatory problems by cement leakage
- Micro movement of the cement against the surface the bone due to inadequate fixation
- Fat embolism
- decrease in blood pressure
- myocardial infarction
- cardiac arrest
- cerebrovascular accident
- pulmonary embolism
- cardiac embolism
- short-term cardiac conduction disorder
- sudden death

### Safety precautions

- Prior to use, thoroughly read these instructions for use and become familiar with the surgical technique.

- Keep the instructions for use accessible to all staff.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrect operating techniques, the limitations of treatment methods or inadequate asepsis.
- Under no circumstances may modular implant components from different suppliers be combined.
- Damage to the weight-bearing structures can give rise to loosening, dislocation and migration, as well as other complications.

Paradigm Spine GmbH | Eisenbahnstrasse 84 | 78573 Wurmlingen | Germany  
Phone +49 (0) 7461-9635990 | Telefax: +49 (0) 7461-963599-20  
info@surgalign.com | www.surgalign.com

LBL 0702 Rev. M 2023-05

- To ensure the earliest possible detection of such catalysts of implant malfunction, the implant must be checked periodically post-operatively using appropriate techniques.
- Only stabilize a segment with either rod couplers or fusion components.
- **Do not use both a rod coupler and a fusion component at the same level!**

- Never reuse an implant. Although the implant may appear undamaged, previous stresses may have created non-visible damage that could result in implant failure.
- Never use implants if the packaging is damaged. An implant with damaged packaging might be damaged itself and thus may not be used.
- Never use implants that are past their expiration date.

## MRI compatibility

The HPS™ implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the HPST™ implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

## Storage, inspection and sterilization

### Storage

The implant is individually packed in protective packaging that is labeled according to its contents.

The implant is sterilized with gamma sterilization (25 kGy minimum).

- Always store the implant in the original protective packaging.
- Do not remove the implant from the packaging until immediately before use.
- Store the implant in a dry and dust-free place (standard hospital environment).

## Disinfection / cleaning

The HPS™ implants are not designed to be disinfected or cleaned by the user.

## Resterilization

The HPS™ implants are not designed to be re-sterilized by the user.

## Technique

The HPS™ implants must be implanted only with the applicable HPST™ instruments. The HPS™ instruments are available from the manufacturer at any time.

## Preoperative

The operating surgeon draws up an operation plan specifying and documenting the following:

- Implant component(s) and their dimensions.
- Proper positioning of the implant components.
- Determination of intra-operative orientation points. The following conditions must be fulfilled prior to application:
- All required implant component(s) are readily available.
- Highly aseptic operating conditions are present.
- All requisite implantation instruments must be available and in working order.

*Never use or process damaged or defective instruments.  
Contact your local representative or dealer for repair or replacement.*

**PARADIGM SPINE**  
*the movement in spine care*

*The use of an instrument for tasks other than those for which they are intended may result in damaged / broken instruments or patient injury.*

- The operating surgeon must be especially trained in spinal surgery, biomechanical principles of the spine and the relevant operating techniques.
- The operating surgeon and operating room team must be thoroughly familiar with the operating technique, as well as the range of implants and instruments to be applied. Complete information on these subjects must be readily available at the workplace.

The operating procedure must be explained to the patient, and the patient's understanding of the following information must be documented:

- The patient is aware of the risks associated with neurosurgery, general surgery, orthopedic surgery and with general anesthesia.
- The patient has been informed about the advantages and disadvantages of the implant procedure and about possible alternative treatments.
- The implant can fail due to excessive load, wear and tear, or infection
- The service life of the implant is determined by body weight and physical activity. The implant must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may be necessary if the implant fails. The patient must have their physician carry out follow-up examinations of the implant at regular intervals.

#### Intra-operative

Prior to use, please read the HPS™ product brochure and become familiar with the surgical technique, the corresponding implants and the instruments.

- Prior to use, verify the integrity of the sterile packaging. Never use implants if the packaging is damaged.
- Prior to use, check the product expiration date. Never use implants that are past their expiration date.

#### Posterior Dynamic Stabilization:

- Observe orientation (marking line) when introducing the implant! (The coupler must be oriented so the laser mark is visible on the upper side of the rod)
- The slots of the rod coupler must remain parallel when tightening the nuts (Note: if the slots of the rod coupler are not parallel, it is necessary to correct the position of the coupler by turning the pedicle screws using the counter torque instrument)
- When using the rod coupler with a rod, damaging the coupler must be avoided, especially during the bending process.

#### Posterior Rigid Stabilization

- Only bend rods with the intended instrument.
- Never bend back pre-bent rods.

#### HPS™ Pedicle screw, fenestrated

If the use of bone cement is required, it has to be introduced via the HPS™ Polyaxial pedicle screw, fenestrated or the HPS™ 2.0 Polyaxial pedicle screws utilizing the HPS™ Cement injection set. The bone cement must have a medium viscosity and is intended only for posterior rigid stabilization.

#### General information regarding technique

- Do not use implants that are damaged with scratches, notches or unintentional bending.
- Only use the instruments intended for the respective step.
- Only use the intended instruments (e.g., bending wrench, bending instruments) to bend the rods.
- Use a torque wrench with a counter torque instrument to avoid bending the coupler and movement of the pedicle screws during tightening.
- The proper tightening torque is reached when the torque wrench indicates the necessary torque.
- Use a distraction and/or compression instrument to obtain sagittal alignment of the implant.

#### NOTE

*Detailed information regarding the use of the HPS™ Hybrid Performance System can be found in the HPS™ product Brochure. Further information regarding the procedure is available in workshops, product training or individual consultation with Paradigm Spine.*

#### Postoperative

- Reiterate preoperative instructions to the patient.
- Ensure that the patient is aware of physical activity restrictions and possible adverse reactions.

#### Revision surgery / implant removal

The HPS™ implants are intended for permanent implantation and are typically not removed. However, removal of the implant can be necessary in the following situations:

- Implant breakage
- Pain due to the implant
- Infection
- Pseudarthrosis
- Allergic reactions

If a revision is performed where rod couplers have been utilized and the surgeon elects to fuse this segment, the rod couplers may be replaced with fusion components. The surgeon shall determine whether new mono and/or polyaxial pedicle screws are used, or if the fusion component should be placed over the existing screws. If it is necessary to replace the rod coupler and/or fusion components, it is advisable to replace them on both sides of the spinal segment. Implants which appear to be intact may have invisible damage.

#### Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with these instructions for use. Subject to technical modifications.

#### For further information

Please contact Paradigm Spine or your medical product representative if further information is needed. The HPS™ System is protected by law worldwide.

#### IMPORTANT NOTE

These instructions for use are complete at the time of publishing.

#### Patent

For patent information see [www.surgalign.com](http://www.surgalign.com)

## Symbols

	Manufacturer
	Date of manufacture
	Use-by date
<b>LOT</b>	Batch code
<b>REF</b>	Catalogue number
<b>QTY</b>	Content of usable units(s)
<b>STERILE R</b>	Sterilized using irradiation
	Do not resterilize
	Do not use if package is damaged
	Keep dry
	Do not reuse
	Consult instructions for use at this website.
<b>MD</b>	Medical Device
	Caution
<b>CE</b> 0123	CE-marking according to directive 93/42/EEC

Paradigm Spine GmbH | Eisenbahnstrasse 84 | 78573 Wurmlingen | Germany  
 Phone +49 (0) 7461-9635990 | Telefax: +49 (0) 7461-963599-20  
 info@surgalign.com | www.surgalign.com

LBL 0702 Rev. M 2023-05

**PARADIGM SPINE**  
*the movement in spine care*

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