coflex® Implant

Package contents

The package contains one coflex® implant.

Description

The coflex® is a dynamic interspinous decompression spacer intended for the lumbar spine in patients with stenosis or mild degenerative instability. The coflex® consists of a single component and is made from titanium alloy (ELI Ti6Al4V). The coflex® has wings designed to provide stability and the "U" shaped body is intended to allow flexion and extension of the lumbar spine.

Material

The coflex® implant is made from wrought titanium 6-aluminum 4-vanadium alloy (ISO 5832-3). It is supplied sterile and is available in numerous sizes.

Intended use

The coflex® implant is intended for permanent implantation between the spinous processes of 1 or 2 lumbar motion segments (L1-L5) and controls segmental motion in cases of lumbar stenosis or mild degenerative instability. The coflex® implant may also be used in up to 2 lumbar motion segments adjacent to fused level(s).

Intended users

The coflex® implant and instruments should only be used by surgeons who are experienced and familiar with the implant components, instruments, procedure, clinical applications, biomechanics, adverse events, and risks associated with coflex®.

Intended patient population

The coflex® implant and instruments are intended to be used in skeletally mature patients.

Indications

Radiographically confirmed moderate to severe stenosis, isolated to 1 or 2 levels, including conditions such as stable grade 1 spondylolisthesis.

Contraindications

The coflex® is contraindicated on patients with:

- Back or leg pain of unknown etiology.
- Axial back pain only, with no leg, buttock, or groin pain
- Cauda equina syndrome defined as neural compression causing neurogenic bowel or bladder dysfunction
- Grade II or greater spondylolisthesis.
- Isthmic spondylolisthesis or spondylolysis (pars fracture).
- Degenerative lumbar scoliosis (Cobb angle of greater than 25°).
- Severe facet hypertrophy that requires extensive bone removal which would cause instability.
- Radiographically compromised vertebral bodies at any lumbar level(s) caused by current or past trauma or tumor (e.g., compression fracture).
- Osteoporosis.
- Prior fusion or decompressive laminectomy at any index lumbar level.
- Active or chronic infection – systemic or local.
- Known allergy to titanium alloys or MR contrasting agents.
- Systemic and metabolic diseases.
- Morbid obesity defined as a body mass index > 40.
- Pregnancy
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Lack of patient cooperation

- Any medical or surgical condition precluding the potential benefit of spinal surgery

Potential Adverse Events

- Allergic reaction for the implant
- Fracture
- Infection
- Paralysis
- Excessive bleeding (no implant placed)
- Restenosis
- Continued pain / lack of effect
- Progression of spinal instability
- Implant breakage
- Bent implant
- Spacer would not deploy, procedure aborted

General Surgery / Procedure related risks

- Wound Infection
- Neurological complications
 - Restenosis/continued pain/lack of effect/Lumbar stenosis symptoms
 - Paralysis
 - Dural tear/damage
- Ectopic ossification
- Osteolysis
- ASD
- Progression of instability
- Cardiovascular and pulmonary complications
 - Excessive bleeding
 - Incision hematoma
- Death
- Urinary tract infection

Device related risks

- Allergic reaction
- Device loosening/migration/misplacement
- Implant failure/breakage/bent implant
- Spinous process fracture
- Revision or reoperation

Intended performance/patient benefits and other device claims

Clinical performance

- coflex® maintained motion at the index level at 5 years. As assessed using the ROM in radiological follow-up images.
- coflex® maintained physiological adjacent segment kinematics at 5 years. As assessed using the ROM in radiological follow-up images.
- coflex® maintained translational motion at the index level at 3 years. As assessed using the ROM in radiological follow-up images.
- coflex® maintained foraminal height at 5 years. As measured on radiological follow-up images
- Treatment with the coflex® resulted in superior foraminal and disc height maintenance (allowing nerves to exit the spinal canal freely) compared to decompression alone.

Clinical performance

- coflex® significantly decreased the degree of back pain (both immediately and long-term).
- coflex® significantly decreased the degree of leg pain (both immediately and long-term).

- coflex® significantly improved physical function.
- coflex® significantly decreased the symptom severity
- More patients were satisfied with the coflex® procedure compared to fusion.
- coflex® patients were more satisfied with their outcomes compared to fusion patients.
- On average, the surgery with coflex® is up to an hour shorter than fusion surgery
- Less time in the operating room, less blood loss, and fewer days in the hospital compared to fusion with pedicle screws
- coflex® patients were able to return home two days earlier compared to fusion patients
- Treatment with the coflex® reduces trauma to your back, decreases the amount of blood loss during surgery, and allows for a shorter hospital or outpatient surgical center stay compared to pedicle screw fusion.
- Treatment with the coflex® resulted in superior composite clinical results (combination of 4 clinically relevant factors: disability improvement, neurological maintenance or improvement, no subsequent interventions, no device-related complications) compared to decompression alone.
- Treatment with the coflex® resulted in superior walking distance improvement (2.4x more compared to decompression alone)
- Treatment with the coflex® resulted in significantly less epidural steroid injections post-surgery (decompression patients had 3.3x more)
- Treatment with the coflex® resulted in faster relief of symptoms and quicker recovery (pain and function measurements) compared to pedicle screw fusion.

Clinical Safety

- For Two level coflex® procedures the revision rate is significantly lower compared to fusion procedure.

Warnings / Precautions

- Prior to use, thoroughly read these instructions for use and become familiar with the surgical technique.
- Keep the instructions for use accessible to all staff.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques. Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrect operating techniques, the limitations of treatment methods or inadequate asepsis.
- Under no circumstances may modular implant components from different suppliers be combined.
- Each patient's record shall document the implant used (name, article number, lot number).
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed about post-surgical behavioral requirements.
- Damage to the weight-bearing structures can give rise to loosening, dislocation and migration, as well as other complications. To ensure the earliest possible detection of such catalysts of implant dysfunction, the implant must be checked periodically post operative using appropriate techniques.

- Fracture of the coflex® implant wing may occur during preparation when the wings are bent open or during insertion when the wings are crimped.

- Never reuse an implant. Although the implant may appear undamaged, previous stresses may have created non-visible damage that could result in implant failure.

- Never use implants if the packaging is damaged. An implant with damaged packaging might be damaged itself and thus may not be used.

- Never use implants that are past their expiration date.

MRI Conditional

Non-clinical testing has demonstrated that the coflex® implant is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5-Tesla (1.5T) or 3.0-Tesla (3.0T)
- Spatial gradient field of up to:
 - 11,230 G/cm (112.3 T/m) for 1.5T systems
 - 5,610 G/cm (56.1 T/m) for 3.0T systems.
- Maximum whole body averaged specific absorption rate (SAR) of:
 - 2.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 1.5T
 - 2.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 3.0T.

3.0T RF heating

In non-clinical testing with body coil excitation, the coflex® Interlaminar Implant produced a temperature rise of less than 3.5 °C at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 3.0T Siemens Trio (MRC20587) MR scanner with SYNGO MR A30 4VA30A software.

1.5T RF heating

In non-clinical testing with body coil excitation, the coflex® Interlaminar Implant produced a temperature rise of less than 3.5 °C at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 1.5T Siemens Espree (MRC30732) MR scanner with SYNGO MR B17 software.

The RF heating behavior does not scale with static field strength. Devices which do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

MR Artifact

In testing using a 3.0T system with spin-echo sequencing, the shape of the image artifact follows the approximate contour of the device and extends radially up to 19 mm from the implant.

Storage, inspection and sterilization

Storage

The implant is individually packed in protective packaging that is labeled according to its contents.

STERILE R The implant is sterilized with gamma sterilization (25 kGy minimum).

- Always store the implant in the original protective packaging.
- Do not remove the implant from the packaging until immediately before use.

- Store the implant in a dry and dust-free place (standard hospital environment).

Disinfection / cleaning

The coflex® implant is not designed to be disinfected or cleaned by the user.

Resterilization

The coflex® implant is not designed to be reesterilized by the user.

Procedure

The coflex® implant must be implanted only with the applicable coflex® trial and instruments. The coflex® trial and instruments are available from the manufacturer at any time.

Preoperative

The operating surgeon draws up an operation plan specifying and documenting the following:

- Implant component(s) and their dimensions.
- Proper position of the implant component(s) between the spinous processes.
- Determination of intraoperative orientation points.

The following conditions must be fulfilled prior to application:

- All required implant component(s) are readily available.
- Highly aseptic operating conditions are present.
- All requisite implantation instruments must be available and in working order, including the coflex® trial

Never use or process damaged or defective instruments. Contact your local representative or dealer for repair or replacement.

The use of an instrument for tasks other than those for which they are intended may result in damaged / broken instruments or patient injury.

- The operating surgeon and operating room team must be thoroughly familiar with the operating technique, as well as the range of implants and instruments to be applied. Complete information on these subjects must be readily available at the workplace.
- The operating surgeon must be especially trained in spinal surgery, biomechanical principles of the spine and the relevant operating techniques.

The operating procedure must be explained to the patient, and the patient's understanding of the following information must be documented:

- The patient is aware of the risks associated with neurosurgery, general surgery, orthopedic surgery and with general anesthesia.
- The patient has been informed about the advantages and disadvantages of the implant procedure and about possible alternative treatments.
- The implant can fail due to excessive load, wear and tear, or infection.
- The service life of the implant is determined by body weight and physical activity. The implant must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may be necessary if the implant fails.
- The patient must have their physician carry out followup examinations of the implant at regular intervals.

Intraoperative

Prior to use, please read the coflex® product brochure and become familiar with the surgical technique, the corresponding implants and the instruments.

- Prior to use, verify the integrity of the sterile packaging. Never use implants if the packaging is damaged.

- Prior to use, check the product expiration date. Never use implants that are past their expiration date.

Perform microsurgical decompression and measure interspinous space utilizing the trial portion of the instrument. Load the coflex® implant on the insertion portion of the instrument and introduce into the prepared site. If desired, the wings of the coflex® implant may be opened with the coflex® bending pliers and closed with the coflex® crimping pliers. Only bend the wings open as much as necessary to safely introduce the implant.

Postoperative

- Reiterate preoperative instructions to the patient.
- Ensure that the patient is aware of physical activity restrictions and possible adverse reactions.

During the first six weeks after surgery the patient should avoid any strenuous physical activity. Examples of strenuous physical activity include lifting more than 5 Kg. The patient should not do sports until the doctor agrees to do so. Sports include swimming, golf, tennis, running, and jogging among others. If the patient is too active too soon after surgery, this may result in implant-related or -unrelated complications postoperatively. This could cause pain. The patient could need further surgery.

The patient should tell the doctor as soon as possible after the surgery if there is pain or swelling in the back or numbness in the legs or buttocks. These symptoms can be a sign that the coflex® implant is not working properly. More surgery may be needed. The patient should inform the doctor about any traumatic event e.g. fall.

Revision surgery / implant removal

The coflex® implant is intended for permanent implantation and is not intended for removal. When device removal is necessary, please follow these instructions:

- Use adequate surgical approach to visualize coflex®.
- The coflex® implant can be removed using standard instruments like chisels, forceps and pliers.
- Care has to be taken to remove bone overgrowing the wings of the device.
- Then the wings can be carefully bent open using a small bone chisel before the device can be taken out using pliers or forceps.

On special request Paradigm Spine offers a revision instrument that is not part of the standard set. Additional procedures may be necessary to stabilize the spine following coflex® implant removal, including implantation of a new coflex® implant or spinal fusion.

Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use. Technical alterations reserved.

For further information

Please contact Paradigm Spine or your authorized representative if further information on this product is needed.

Patent




For patent information see www.paradigmspine.com

Summary of safety and clinical performance

For information see www.paradigmspine.com

Symbols

	Medical Device
	Manufacturer
	Date of manufacture
	Use-by date
	Lot number
	Catalog number
	Content of usable units(s)
	Sterilized using irradiation
	Do not resterilize
	Do not use if package is damaged
	Keep dry
	Do not reuse
	Consult instructions for use at this website.
	Caution
	MR conditional
	CE-marking according to Regulation (EU) 2017/745
	CE-marking according to Regulation (EU) 2017/745
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Unique Device Identification

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