

Butrex[®]

Lumbar Buttress Plating System



- › Resilient locking arm technology
- › Serrated bottom surface for greater grip
- › Monoblock construction


Butrex[®]

Lumbar Buttress Plating System

Item Number	Description
X023-0001	22 mm Butrex Plate
X023-0002	26 mm Butrex Plate
X023-0011	4.3 mm Fixed Screw, 20 mm
X023-0012	4.3 mm Fixed Screw, 22 mm
X023-0013	4.3 mm Fixed Screw, 24 mm
X023-0021	4.5 mm Fixed Screw, 20 mm
X023-0022	4.5 mm Fixed Screw, 22 mm
X023-0023	4.5 mm Fixed Screw, 24 mm

Monoblock Construction
Eliminates post-op disassembly due to fewer parts

Pre-contoured Plate
Conforms to local anatomy



Graft Viewing Window
Provides ideal graft/endplate interface visualization

Low Profile/Volume Plate
Minimizes tissue disruption

**Self-Drilling/
Self-Tapping Screws**
Ease of insertion

Serrated Surface
Reduces risk of interoperative plate migration

Two Point Fixation
Higher cantilever
expulsion resistance



888.886.9354
cs@xtantmedical.com
xtantmedical.com
 664 Cruiser Ln,
 Belgrade, MT 59714

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.
 LIMITED WARRANTY and DISCLAIMER: Xtant Medical products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed. WARNING: In the USA, this product has labeling limitations. See package insert for complete information. CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

Butrex[®] is a trademark of Xtant Medical.
 © Xtant Medical. All Rights Reserved. FM-C-MRK-53 Rev C 11/21