


STREAMLINED VERSATILITY




Streamline® OCT
Occipito-Cervico-Thoracic System




Versatile Occipital Plate

- Modular system provides cephalad-caudal, sagittal, coronal and medial-lateral adjustability, reducing the need for extreme rod bends




Multiple Screw Options

- Polyaxial, smooth shank and high angle
 - Standard diameters: 3.5, 4.0, 4.5mm
 - Standard lengths: 36, 42, 48mm
- Low-profile friction fit screw heads
- High angle screw allows for 360° positioning and up to 60° angulation



Wide Variety of Connectors

- Rod-to-rod connectors, 3.2-3.2mm, 3.2-5.5mm
- Parallel, straight and angled lateral offset connectors with closed and open rod attachments



Crosslink Connectors

- Rod-to-rod and head-to-head connectors with 8mm of medial/lateral adjustability

Comprehensive Rod Offering

- Low profile 3.2mm diameter CoCr
- Straight, pre-bent and transition

Contact your RTI Surgical representative for complete labeling and product information.



RTI Surgical, Inc.

® indicates U.S. trademark registration. All trademarks and/or images are the property of their respective owners or holders.

Distributed by:
RTI Surgical, Inc.
375 River Park Circle
Marquette, MI 49855 USA
t: 800.557.9909
www.rtisurgical.com