


Cortera

Posterior Fixation System




Thoracolumbar Pedicle Screw System

Cortera

Posterior Fixation System

The Cortera Spinal Fixation System is a comprehensive posterior thoracolumbar fixation solution. Featuring innovative implants and multi-functional instrumentation, Cortera provides surgeons with a safe and effective solution that can improve surgical work flow and deliver value when navigating complex procedures.

- Instrumentation designed with features to easily handle straightforward constructs and simplify complex cases
- Dual lead and Dual-to-Quad Lead Screw Thread Options
- 5.5/6.0mm, Ti/CoCr Rod System
- Helical Flange Design
- Universal S27 Hexalobe Drive



☎ 888.886.9354
✉ cs@xtantmedical.com
🌐 xtantmedical.com
📍 664 Cruiser Ln,
Belgrade, MT 59714

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.
LIMITED WARRANTY and DISCLAIMER: Xtant Medical products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed. WARNING: In the USA, this product has labeling limitations. See package insert for complete information. CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

Cortera is a trademark of Xtant Medical
© Xtant Medical. All Rights Reserved. MRK-000142 (A) 9/2024